



Discovery™

Elbow System

POLYETHYLENE EXCHANGE

SURGICAL TECHNIQUE

 **Lima Corporate**
Orthopaedic  motion

limacorporate.com



DISCOVERY™ ELBOW POLYETHYLENE EXCHANGE

One Surgeon. One Patient.®

Over 1 million times per year, Biomet helps one surgeon provide personalized care to one patient.

The science and art of medical care is to provide the right solution for each individual patient. This requires clinical mastery, a human connection between the surgeon and the patient, and the right tools for each situation.



At Biomet, we strive to view our work through the eyes of one surgeon and one patient. We treat every solution we provide as if it's meant for a family member.

Our approach to innovation creates real solutions that assist each surgeon in the delivery of durable personalized care to each patient, whether that solution requires a minimally invasive surgical technique, advanced biomaterials or a patient-matched implant.

When one surgeon connects with one patient to provide personalized care, the promise of medicine is fulfilled.

DISCOVERY™ ELBOW POLYETHYLENE EXCHANGE

Preoperative Planning\Bearing Removal



▼ PREOPERATIVE PLANNING

Place two Ulna Bearing Revision kits (part number 114800) in a freezer for a minimum of three hours. Freezing causes the bearing to constrict, making it easier to insert into the ulna ring. The temperature should be between -13° F and 14° F (-25° C and -10° C). A lower freezer temperature will increase the handling time of the bearing. Do not remove the bearing from the freezer until ready for assembly, as it will begin to expand immediately and reach full expansion within two minutes of removal. The second bearing kit should remain in the freezer as reserve.

Caution: Read through the entire bearing insertion kit technique prior to removing bearing from the freezer.

▼ BEARING REMOVAL

Insert the threaded end of the bearing removal tool through the center hole of the ulna component. Push the T-handle toward the ulna component as far as possible. If necessary, turn the T-handle to allow the threaded shaft to pass through the bearing. The ledge of the ulna component should fit into the recess on the body of the bearing removal tool (Figure 1). With the jagged end toward the polyethylene, tighten the end cap onto the bearing removal tool (Figure 2).

This brochure is presented to demonstrate the surgical technique of Hill Hastings II, M.D. Biomet, as the manufacturer of this device, does not practice medicine and does not recommend this device or technique. Each surgeon is responsible for determining the appropriate device and technique to utilize on each individual patient.

DISCOVERY™ ELBOW POLYETHYLENE EXCHANGE

Bearing Removal \New Bearing Insertion




figure 3



figure 4

▼ BEARING REMOVAL (CONT.)

While holding the body of the removal tool, rotate the T-handle clockwise. The polyethylene will be pulled from within the ulna component onto the threaded shaft of the removal tool. Continue rotating until the polyethylene is removed from the ulna stem and the locking pin falls free (Figure 3). Discard the pin and polyethylene; irrigate and remove any small polyethylene particles

▼ NEW BEARING INSERTION

Remove one ulna bearing revision kit from the freezer. The widest portion of the bearing should face toward the widest portion of the ulna ring (Figure 4). Locate the four notches on the outer edge of the bearing and align the cylindrical notch (pin groove) posteriorly. Align the three shallow notches in the bearing with the

DISCOVERY™ ELBOW POLYETHYLENE EXCHANGE

New Bearing Insertion\Optional Bearing Insertion Method



Figure 5



Figure 6

three tabs inside the ulna stem ring (Figure 5). Insert the bearing into the middle of the stem ring by pushing it until the bearing freely spins/rotates. Rotate the bearing until the cylindrical notch (pin groove) aligns with the pin groove of the ulna stem (Figure 6).



Figure 7



Figure 8

▼ OPTIONAL BEARING INSERTION METHOD

The bearing rotation tool may be used to rotate the bearing if it does not spin freely when inserted into the ulna ring.

To use, insert the tool from the lateral side of the bearing with the correct anatomic engraving facing outward, allowing the long metal tab to slide into the groove on the bearing reserved for the locking pin (Figure 7).

Place the Discovery™ screw driver through the drive hole at the end of the tool (this will aid in rotation). Rotate the bearing until the cylindrical notch (pin groove) aligns with the pin groove of the ulna stem (Figure 8). If done correctly, the long metal tab will stop against the widest portion of the ulna ring.

DISCOVERY™ ELBOW POLYETHYLENE EXCHANGE

Pin Insertion



Figure 9



Figure 10




Figure 11




Figure 12

▼ PIN INSERTION

Place the ulna pin inserter guide flush against the ulna stem, ensuring proper alignment of the guide tube with the pin hole in the ulna stem. Insert the pin with the chamfered end first (Figure 9 and 10) through the guide tube in the pin inserter.

Note: This can be accomplished from either the lateral or medial side.







Thread the pin inserter plunger over the guide tube and thread by hand until resistance is met (Figure 11). Provide downward force on the pin inserter guide to keep proper alignment.

Place the ratchet handle onto the pin inserter plunger and thread the plunger in until the pin is fully seated. The pin is seated when it is flush against the ulna stem (Figure 12).

Caution: Do **NOT** overtighten or use with power.

DISCOVERY™ ELBOW POLYETHYLENE EXCHANGE

Ordering Information

PRODUCT	PART NUMBER	DESCRIPTION	SIZE
	114800*	Discovery™ Ulna Bearing Revision Kit	-
	405908**	Driver Ratchet Handle	-
	414922**	Screwdriver Handle	2.0/2.7mm
	414923**	X-lock Standard Blade (Screwdriver Shaft)	2.4 mm
	414926**	T10 Hexalobular Drive Blade	-
	414950**	Discovery™ Bearing Removal Tool	-
	DSC-419601** DSC-419602**	Discovery™ Bearing Ulna Pin Inserter Guide Discovery™ Bearing Ulna Pin Inserter Plunger	-
	DSC-419600**	Discovery™ Bearing Rotation Tool	-

*Contains one polyethylene bearing and locking pin

**Available from loaners

Limacorporate spa
Via Nazionale, 52
33038 Villanova di San Daniele
Udine - Italy
Tel.: +39 0432 945511
Fax: +39 0432 945512
E-mail: info@limacorporate.com
www.limacorporate.com

Lima Implantes slu
Fonsanta, 46 5ª planta
08970 Sant Joan Despi
Barcelona - Spain
Tel.: +34 93 480 8505
Fax: +34 93 419 6527

Lima France sas
Les Espaces de la Sainte Baume
Parc d'Activité de Gémenos - Bât.A5
30 Avenue du Château de Jouques
13420 Gémenos - France
Tel.: +33 (0) 4 42 01 63 12
Fax: +33 (0) 4 42 04 17 25
E-mail: info@limafrance.com

Lima O.I. doo
Ante Kovacic 3
10000 Zagreb - Croatia
Tel.: +385 (0) 1 2361 740
Fax: +385 (0) 1 2361 745
E-mail: lima-oi@lima-oi.hr

Lima Switzerland sa
Birkenstrasse, 49
CH-6343 Rotkreuz - Zug
Switzerland
Tel.: +41 (0) 41 747 06 60
Fax: +41 (0) 41 747 06 69
E-mail: info@lima-switzerland.ch

Lima Japan kk
Shinjuku Center Building, 29th floor
1-25-1, Nishi-shinjuku, Shinjuku,
Tokyo - Japan
163-0629
Tel.: +81 3 5322 1115
Fax: +81 3 5322 1175

Lima CZ sro
Do Zahrádek I., 157/5
155 21 Praha 5 - Zličín
Czech Republic
Tel.: +420 222 720 011
Fax: +420 222 723 568
E-mail: info@limacz.cz

Lima Deutschland GmbH
Kapstadtring 10
22297 Hamburg - Germany
Tel.: +49 40 6378 4640
Fax: +49 40 6378 4649
E-mail: info@lima-deutschland.com

Lima Austria GmbH
Seestadtstrasse 27 / Top 6-7
1220 Wien - Austria
Tel.: +43 (1) 2712 469
Fax: +43 (1) 2712 469 100
E-mail: office@lima-austria.at

Lima SK s.r.o.
Zvolenská cesta 14
97405 Banská Bystrica - Slovakia
Tel.: +421 484 161 133
Fax: +421 484 161 138
E-mail: info@lima-sk.sk

Lima Netherlands
Havenstraat 30
3115 HD Schiedam
The Netherlands
Tel.: +31 (0) 10 246 26 60
Fax: +31 (0) 10 246 26 61
E-mail: info@limanederland.nl
www.limanederland.nl

Lima Implantes Portugal S.U. Lda
Rua Olavo D'Eça Leal Nº6 Loja-1
1600-306 Lisboa - Portugal
Tel.: +35 121 727 233 7

Lima Orthopaedics Australia Pty Ltd
Unit 1, 40 Ricketts Rd
Mt Waverley 3149
Victoria Australia
Tel.: +61 (03) 9550 0200
Fax: +61 (03) 9543 4003
www.limaortho.com.au

Lima Orthopaedics New Zealand Ltd
Zone 23, Unit 102, Edwin Street, Mt Eden
Auckland, 1024
New Zealand
Tel.: +64 (09) 531 5522
Fax: +64 (09) 522 3380

Lima Orthopaedics UK Limited
Unit 2, Campus 5
Third Avenue
Letchworth Garden City
Herts. SG6 2JF
United Kingdom
Tel.: 0844 332 0661
Fax: 0844 332 0662

Lima USA Inc.
2001 NE Green Oaks Blvd., Suite 100
Arlington, TX 76006
Tel.: +1 817-385-0777
Fax: +1 817-385-0377

Lima Sweden AB
Företagsallén 14 B
SE-184 40 Akersberga
Sweden
Tel.: +46 8 544 103 87
Fax: +46 8 540 862 68
www.linksSweden.se

Lima Italy
Centro Direzionale Milanofiori
Strada 1 - Palazzo F1
20090 Assago - Milano
Italy
Tel.: +39 02 57791301

Lima Korea Co. Ltd.
Zero Bldg., 11FL,
14 Teheran-ro 84-gil, Gangnam Gu,
Seoul, 135-845, Korea
Tel.: +82-2-538-4212
Fax: +82-2-538-0706

Lima do Brasil Ltda
Av. Sagitário, 138 - Sala 2707
Edifício Torre City - Condomínio Alpha Square
06473-073 - Barueri - SP - Brasil
Tel.: +55 (11) 2664-0620
Fax: +55 (11) 2664-0621
E-mail: info@limadobrasil.com.br
www.limadobrasil.com.br

Lima SM spa
Strada Borrana 38
47899 Serravalle, Republic of San Marino
Tel.: +378 0549 961911
Fax: +378 0549 961912
E-mail: info@lima-sm.com

Discovery™ Trademark is licensed to Limacorporate S.P.A. for use only in the European Economic Area (E.E.A.) and Switzerland.

For product information, including indications, contraindications, warnings, precautions and potential adverse effects, see the package insert or Biomet's website.

limacorporate.com

This publication is not intended for distribution in the U.S.

The DISCOVERY™ Elbow System is manufactured by Biomet.

The DISCOVERY™ Elbow System is distributed in the European Economic Area (E.E.A.) and Switzerland by Limacorporate spa

B.1402.23.000.1

111500

